

Product Datasheet

PEX3 (IN) Peptide (orb1231633)

Description	PEX3 (IN) Peptide
Tested Applications	Blocking
Target	PEX3
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store PEX3 (IN) Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	PEX3 (IN) Peptide is used for blocking the activity of the PEX3 (IN) Antibody
Uniprot ID	P56589
NCBI	P56589
Expiration Date	6 months from date of receipt.