

Product Datasheet

GAPDH [12D3H9] Peptide (orb1231655)

Description	GAPDH [12D3H9] Peptide
Tested Applications	Blocking
Target	GAPDH
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store GAPDH peptide at -20°C, stable for one year.
Note	For research use only
Application notes	GAPDH peptide is used for blocking the activity of GAPDH antibody.
Uniprot ID	P04406
NCBI	NP_002037
Expiration Date	6 months from date of receipt.