

Product Datasheet

SARS-CoV-2 (COVID-19) Alpha Variant (B.1.1.7, UK) Spike P681H Peptide (orb1231664)

Description	SARS-CoV-2 (COVID-19) Alpha Variant (B.1.1.7, UK) Spike P681H Peptide
Tested Applications	Blocking
Target	S
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store SARS-CoV-2 (COVID-19) Alpha Variant (B.1.1.7, UK) Spike P681H at -20°C, stable for one year.
Note	For research use only
Application notes	SARS-CoV-2 (COVID-19) Alpha Variant (B.1.1.7, UK) Spike P681H is used for blocking the activity of the SARS-CoV-2 (COVID-19) Alpha Variant (B.1.1.7, UK) Spike P681H Antibody.
NCBI	QHD43416
Expiration Date	6 months from date of receipt.