

Product Datasheet

SARS-CoV-2 (COVID-19) NSP13 (Helicase) Peptide (orb1231679)

Description	SARS-CoV-2 (COVID-19) NSP13 (Helicase) Peptide
Tested Applications	Blocking
Target	ORF1ab
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store SARS-CoV-2 (COVID-19) NSP13 (Helicase) Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	SARS-CoV-2 (COVID-19) NSP13 (Helicase) Peptide is used for blocking the activity of the SARS-CoV-2 (COVID-19) NSP13 (Helicase) Antibody.
Uniprot ID	P0DTC1
NCBI	YP_009725308.1
Expiration Date	6 months from date of receipt.