


Product Datasheet

UBD Peptide (orb1231779)

Description	UBD Peptide
Tested Applications	Blocking
Target	UBD
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store UBD peptide at -20°C, stable for one year.
Note	For research use only
Application notes	UBD peptide is used for blocking the activity of UBD antibody.
Uniprot ID	O15205
NCBI	NP_006389
Expiration Date	6 months from date of receipt.



Western blot analysis of UBD in human co...