

Product Datasheet

SNCA Peptide (orb1231878)

Description	SNCA Peptide
Tested Applications	Blocking
Target	SNCA
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store SNCA Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of SNCA Antibody (catalog number 7911). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P37840
NCBI	NP_000336
Expiration Date	6 months from date of receipt.