

Product Datasheet

TAS2R4 antibody (orb12319)

Description	Goat polyclonal antibody to TAS2R4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	TAS2R4 (aa157-171)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from
Protein Sequence	PELVTTRNNTSFNIS
Clonality	Polyclonal
MW	33.8
Entrez	50832
Dilution Range	ELISA: 1:128000
Expiration Date	12 months from date of receipt.