

Product Datasheet

PPARGC1A Peptide (orb1231913)

Description	PPARGC1A Peptide
Tested Applications	Blocking
Target	PPARGC1A
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store PPARGC1A Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of PPARGC1A Antibody (catalog number 7793). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q9UBK2
NCBI	NP_037393
Expiration Date	6 months from date of receipt.