

Product Datasheet

Sin Nombre Virus Glycoprotein 2 Peptide (orb1231957)

Description	Sin Nombre Virus Glycoprotein 2 Peptide
Tested Applications	Blocking
Target	SNVsMgp1
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store Sin Nombre Glycoprotein 2 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of Sin Nombre Virus Glycoprotein 2 antibody (catalog number 7683). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
NCBI	NP_941974
Expiration Date	6 months from date of receipt.