

## Product Datasheet

### Andes Virus Glycoprotein Peptide (orb1231960)

|                            |                                                                                                                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Andes Virus Glycoprotein Peptide                                                                                                                                                                                                         |
| <b>Tested Applications</b> | Blocking                                                                                                                                                                                                                                 |
| <b>Target</b>              | ANDVsMgp1                                                                                                                                                                                                                                |
| <b>Preservatives</b>       | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide                                                                   |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                                   |
| <b>Concentration</b>       | 200 ug/mL                                                                                                                                                                                                                                |
| <b>Storage</b>             | Store Andes Virus Glycoprotein peptide at -20°C, stable for one year.                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                    |
| <b>Application notes</b>   | The peptide is used for blocking the activity of Andes Virus Glycoprotein antibody (catalog number 7677). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C. |
| <b>NCBI</b>                | <a href="#">AAO86638</a>                                                                                                                                                                                                                 |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                           |