

Product Datasheet

SPINSTER Peptide (orb1231961)

Description	SPINSTER Peptide
Tested Applications	Blocking
Target	SPNS1
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store SPINSTER peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of SPINSTER antibody (catalog number 7651). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q9H2V7
NCBI	NP_001135920
Expiration Date	6 months from date of receipt.