

Product Datasheet

CCR7 Peptide (orb1231962)

Description	CCR7 Peptide
Tested Applications	Blocking
Target	CCR7
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store CCR7 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of CCR7 antibody (catalog number 7649). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P32248
NCBI	NP_001829
Expiration Date	6 months from date of receipt.