

Product Datasheet

CXCR3 Peptide (orb1231963)

Description	CXCR3 Peptide
Tested Applications	Blocking
Target	CXCR3
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store CXCR3 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of CXCR3 antibody (catalog number 7647). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P49682
NCBI	NP_001136269
Expiration Date	6 months from date of receipt.