


Product Datasheet

ITCH Peptide (orb1231964)

Description	ITCH Peptide
Tested Applications	Blocking
Target	ITCH
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store ITCH peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of ITCH antibody (catalog number 7645). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q96J02
NCBI	NP_001244066
Expiration Date	6 months from date of receipt.



Western blot analysis of ITCH in 3T3 cel...