

## Product Datasheet

### FLI1 Peptide (orb1232025)

|                            |                                                                                                                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | FLI1 Peptide                                                                                                                                                                                                         |
| <b>Tested Applications</b> | Blocking                                                                                                                                                                                                             |
| <b>Target</b>              | FLI1                                                                                                                                                                                                                 |
| <b>Preservatives</b>       | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide                                               |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                               |
| <b>Concentration</b>       | 200 ug/mL                                                                                                                                                                                                            |
| <b>Storage</b>             | Store FLI1 peptide at -20°C, stable for one year.                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                |
| <b>Application notes</b>   | The peptide is used for blocking the activity of FLI1 antibody (catalog number 7449). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C. |
| <b>Uniprot ID</b>          | <b>P20930</b>                                                                                                                                                                                                        |
| <b>NCBI</b>                | <b>NP_002008</b>                                                                                                                                                                                                     |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                       |