

Product Datasheet

APO-E Peptide (orb1232061)

Description	APO-E Peptide
Tested Applications	Blocking
Target	APOE
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store APO-E peptide at -20°C, stable for one year.
Note	For research use only
Application notes	APO-E peptide is used for blocking the activity of APO-E antibody (catalog number 7333). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P02649
NCBI	NP_000032
Expiration Date	6 months from date of receipt.