

Product Datasheet

MTERFD3 Peptide (orb1232075)

Description	MTERFD3 Peptide
Tested Applications	Blocking
Target	MTERFD3
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store MTERFD3 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	MTERFD3 peptide is used for blocking the activity of MTERFD3 antibody (catalog number 7285). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q49AM1
NCBI	NP_079474
Expiration Date	6 months from date of receipt.