

Product Datasheet

MTERFD2 Peptide (orb1232076)

Description	MTERFD2 Peptide
Tested Applications	Blocking
Target	MTERFD2
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store MTERFD2 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	MTERFD2 peptide is used for blocking the activity of MTERFD2 antibody (catalog number 7283). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q7Z6M4
NCBI	NP_872307
Expiration Date	6 months from date of receipt.