



## Product Datasheet

### Fubp1 antibody (orb12321)

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal to Fubp1.                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Reactivity</b>          | Bovine, Canine, Human, Mouse, Porcine, Rat                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target</b>              | Fubp1 (mouse, aa160-174)                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                                                                                                 |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Application notes</b>   | WB: Approx 75kDa band observed in lysates of cell line Jurkat (calculated MW of 67.4kDa according to NP_476513.2). Recommended concentration: 0.03-0.1µg/ml.. IHC: In paraffin embedded Human Tonsil shows staining in a selection of cells mostly outside the germinal centre. Recommended concentration: 4-8µg/ml.. Notes: Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371. |
| <b>Protein Sequence</b>    | DQIVEKGRPAPGFHH                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>MW</b>                  | 67.4                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Entrez</b>              | <b>8880</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Dilution Range</b>      | ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 4-8 µg/ml                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                                             |



Western blot analysis of Jurkat lysate u...



Immunohistochemical staining of Human To...