


Product Datasheet

Syt5 antibody (orb12322)

Description	Goat polyclonal to Syt5.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Synaptotagmin-5 (aa 80-93)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human Brain (Olfactory bulb) lysates (calculated MW of 43.1kDa according to NP_058604.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	SYIDKVQPEIEELD
Clonality	Polyclonal
MW	43.1
Entrez	6861
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Olfactory bulb

...