


Product Datasheet

SYT9 antibody (orb12323)

Description	Goat polyclonal to SYT9.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Synaptotagmin-9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Mouse Liver, Pancreas and Thymus lysates (calculated MW of 56.3kDa according to NP_068689.2). Recommended concentration: 0.5-2µg/ml.
Protein Sequence	EHDSCQDFIYHLRD
Clonality	Polyclonal
MW	56.2
Entrez	143425
Dilution Range	ELISA: 1:64000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Liver lys...