


Product Datasheet

DPAGT1 Peptide (orb1232338)

Description	DPAGT1 Peptide
Tested Applications	Blocking
Target	DPAGT1
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store DPAGT1 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	DPAGT1 peptide is used for blocking the activity of DPAGT1 antibody.
Uniprot ID	Q9H3H5
NCBI	NP_001373
Expiration Date	6 months from date of receipt.



Western blot analysis of DPAGT1 in mouse...