

Product Datasheet

MYBPC2 Peptide (orb1232624)

Description	MYBPC2 Peptide
Tested Applications	Blocking
Target	MYBPC2
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store MYBPC2 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	MYBPC2 peptide is used for blocking the activity of MYBPC2 antibody.
Uniprot ID	Q14324
NCBI	EAW71866
Expiration Date	6 months from date of receipt.