


Product Datasheet

ABHD14B antibody (orb12328)

Description	Goat polyclonal to ABHD14B.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ABHD14B (aa188-200)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical protein: NP_116139.1, NP_001139786.1. This antibody is expected to NOT cross react with ABHD14A.
Protein Sequence	HPCYLDKPEEWHT
Clonality	Polyclonal
MW	22.3
Entrez	84836
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human, Mouse, R...