

## Product Datasheet

### Stk39 Peptide (orb1232919)

|                            |                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Stk39 Peptide                                                                                                                                                          |
| <b>Tested Applications</b> | Blocking                                                                                                                                                               |
| <b>Target</b>              | STK39                                                                                                                                                                  |
| <b>Preservatives</b>       | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                 |
| <b>Concentration</b>       | 200 ug/mL                                                                                                                                                              |
| <b>Storage</b>             | Store Stk39 peptide at -20°C, stable for one year.                                                                                                                     |
| <b>Note</b>                | For research use only                                                                                                                                                  |
| <b>Application notes</b>   | Stk39 peptide is used for blocking the activity of Stk39 antibody.                                                                                                     |
| <b>Uniprot ID</b>          | <b>Q9UEW8</b>                                                                                                                                                          |
| <b>NCBI</b>                | <b>NP_037365</b>                                                                                                                                                       |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                         |