


Product Datasheet

Slc22A17 Peptide (orb1233427)

Description	Slc22A17 Peptide
Tested Applications	Blocking
Target	SLC22A17
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store Slc22A17 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	Slc22A17 peptide is used for blocking the activity of Slc22A17 antibody.
Uniprot ID	Q8WUG5
NCBI	Q8WUG5
Expiration Date	6 months from date of receipt.



Western blot analysis of Slc22A17 in SK-...