


Product Datasheet

Sox9b antibody (orb12346)

Description	Goat polyclonal to Sox9b.
Species/Host	Goat
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Sox9b (zebrafish aa5-18)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	No cross reactivity with Sox9a
Protein Sequence	QRGLKMSVSGAPSP
Clonality	Polyclonal
MW	45
Dilution Range	ELISA: 1:32000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Prostate tissue showing Sox9b expression.