

Product Datasheet

Cuff antibody (orb12348)

Description	Goat polyclonal to Cuff.
Species/Host	Goat
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA
Target	cutoff / cuff (Drosophila melanogaster)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from
Protein Sequence	RRDENAHVHNIEKT
Clonality	Polyclonal
MW	44.8
Dilution Range	ELISA: 1:128000
Expiration Date	12 months from date of receipt.