

Product Datasheet

PJA1 antibody (orb12350)

Description	Goat polyclonal antibody to PJA1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Praja-1 / RNF70 (aa285-298)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 30kDa band in Human Brain and Mouse Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 50.4kD
Protein Sequence	HSTSRWRDTANDNE
Clonality	Polyclonal
MW	71.0; 64.7; 50.4
Entrez	64219
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.