

Product Datasheet

SDCCAG8 antibody (orb12352)

Description	Goat polyclonal to SDCCAG8.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SDCCAG8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in lysates of cell line HEK293 (calculated MW of 82.7kDa according to NP_006633.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.
Protein Sequence	TINDQSQYIHHLEAE
Clonality	Polyclonal
MW	82.7
Entrez	10806
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.

