


Product Datasheet

Trp63 antibody (orb12357)

Description	Goat polyclonal to Trp63.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Trp63 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize reported isoforms d, e, f, g, h. NP_035771.1;NP_001120736.1;NP_001120734.1;NP_001120735.1;NP_001120737.1.
Protein Sequence	ENNAQTQFSEPQYC
Clonality	Polyclonal
MW	65.8; 65.3
Entrez	8626
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 lysate using Trp63 antibody (orb12357).