



Product Datasheet

USP6 antibody (orb12361)


Description	Goat polyclonal to USP6.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	USP6 (aa142-155)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Testis, Prostate and Placenta lysates gave no specific signal but low background at antibody concentration up to 1µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration:
Protein Sequence	HHIDL DVRTTLRNH
Clonality	Polyclonal
MW	159
Entrez	9098
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A431 cell...



Immunofluorescence analysis of U2OS cell...



Flow cytometric analysis of A431 cells u...