

Product Datasheet

PLEKHA5 antibody (orb12362)

Description	Goat polyclonal to PLEKHA5.
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	PLEKHA5 / PEPP2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoform 1 NP_061885.2 and isoform 2 NP_001137293.2 only.
Protein Sequence	RENDVKPDHETP
Clonality	Polyclonal
MW	127; 134; 8.72
Entrez	54477
Dilution Range	ELISA: 1:2000, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.