

Product Datasheet

Lipe Peptide (orb1236458)

Description	Lipe Peptide
Tested Applications	Blocking
Target	LIPE
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store Lipe peptide at -20°C, stable for one year.
Note	For research use only
Application notes	Lipe peptide is used for blocking activity of lipe antibody.
Uniprot ID	Q05469
NCBI	ABA03168
Expiration Date	6 months from date of receipt.