

Product Datasheet

Slug Peptide (orb1236460)

Description	Slug Peptide
Tested Applications	Blocking
Target	SNAI2
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store Slug peptide at -20°C, stable for one year.
Note	For research use only
Application notes	Slug peptide is used for blocking activity of Slug antibody.
Uniprot ID	043623
NCBI	043623
Expiration Date	6 months from date of receipt.