

## Product Datasheet

### TOLLIP Peptide (orb1236529)

|                            |                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | TOLLIP Peptide                                                                                                                                                         |
| <b>Tested Applications</b> | Blocking                                                                                                                                                               |
| <b>Target</b>              | TOLLIP                                                                                                                                                                 |
| <b>Preservatives</b>       | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                 |
| <b>Concentration</b>       | 200 ug/mL                                                                                                                                                              |
| <b>Storage</b>             | Store TOLLIP peptide at -20°C, stable for one year.                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                  |
| <b>Application notes</b>   | TOLLIP peptide is used for blocking the activity of TOLLIP antibody.                                                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q9H0E2</a>                                                                                                                                                 |
| <b>NCBI</b>                | <a href="#">AAH18272</a>                                                                                                                                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                         |