

Product Datasheet

TOLLIP Peptide (orb1236530)

Description	TOLLIP Peptide
Tested Applications	Blocking
Target	TOLLIP
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store TOLLIP peptide at -20°C, stable for one year.
Note	For research use only
Application notes	TOLLIP peptide is used for blocking the activity of TOLLIP antibody.
Uniprot ID	Q9H0E2
NCBI	AAH18272
Expiration Date	6 months from date of receipt.