

Product Datasheet

Influenza A H5N1 (Avian) Neuraminidase, C-Terminal Peptide (orb1236901)

Description	Influenza A H5N1 (Avian) Neuraminidase, C-Terminal Peptide
Tested Applications	Blocking
Preservatives	PBS, pH 7.2 containing 0.1% rAlbumin.
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Upon receipt, store (up to one year) at 2-8°C.
Note	For research use only
Application notes	Suitable for use in blocking the activity of the anti-avian influenza A (H5N1) neuraminidase antibody . It usually blocks the antibody activity completely in ELISA assays by incubating the peptide with equal volume of antibody for 30 minutes at 37°C. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
NCBI	CAC95053
Expiration Date	6 months from date of receipt.