

Product Datasheet

Control Normal Rat IgG2a [KLH/G2a-1-1] (Biotin) Antibody (orb1237087)

Description	Control Normal Rat IgG2a [KLH/G2a-1-1] (Biotin) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, IP
Preservatives	Supplied in PBS/sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Normal Rat IgG2a (Biotin) can be used as control in flow cytometry and immunohistochemistry assays.
Isotype	IgG2a
Dilution Range	Normal Rat IgG2a (Biotin) can be used as control in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.