

Product Datasheet

Control Normal Rat IgG1 [KLH/G1-2-2] (Biotin) Antibody (orb1237097)

Description	Control Normal Rat IgG1 [KLH/G1-2-2] (Biotin) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC
Preservatives	Supplied in PBS/sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Normal Rat IgG1 (Biotin) can be used as control in flow cytometry and immunocytochemistry assays.
Isotype	IgG1
Dilution Range	Normal Rat IgG1 (Biotin) can be used as control in flow cytometry and immunocytochemistry assays.
Expiration Date	12 months from date of receipt.