


Product Datasheet

Natalizumab (ITGA4/CD49d) - Research Grade Biosimilar Antibody (orb1237512)

Description	Natalizumab (ITGA4/CD49d) - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / ITGA4/CD49d [Homo sapiens]
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Isotype	IgG4-kappa
Clonality	Recombinant
Source	CHO cells
Uniprot ID	P13612
Expiration Date	12 months from date of receipt.



SDS-
PAGE
Image