

Product Datasheet

Indatuximab Ravtansine (SDC1/CD138) - Research Grade Biosimilar Antibody (orb1237620)

Description Indatuximab Ravtansine (SDC1/CD138)
- Research Grade Biosimilar Antibody

Conjugation Unconjugated

Immunogen Chimeric / SDC1/CD138 [Homo sapiens]

Preservatives PBS buffer pH 7.5

Concentration batch dependent

Storage Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.


Note For research use only

Isotype IgG4-kappa

Clonality Recombinant

Source CHO cells

Expiration Date 12 months from date of receipt.



SDS-
PAGE
Image