

Product Datasheet

C5 Biosimilar Antibody (orb1237911)

Description	C5 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Chimeric / C5
Target	C5
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG4,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2250440-41-4
Uniprot ID	P01031
NCBI	P01031
Expiration Date	12 months from date of receipt.