



Product Datasheet

TBP antibody (orb12380)


Description	Goat polyclonal antibody to TBP
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TBP /Transcription factor IID (aa39-50)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in nuclear, but not in cytosolic lysates of cell line HeLa (calculated MW of 37.7kDa according to NP_003185.1) and 37kDa observed in lysates of NIH3T3 and of Mouse and Rat Testis lysates (calculated MW of 34.7kDa according to Mo
Protein Sequence	TPQPIQNTNSLS
Clonality	Polyclonal
MW	37.7; 35.7
Entrez	6908
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of cel...



Western blot analysis of staining of HeL...



Western blot analysis of Pig Testis lysa...