

Product Datasheet

MS4A1 Biosimilar Antibody (orb1238038)

Description	MS4A1 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / MS4A1 (CD20, MS4A-1) [Homo sapiens]
Target	MS4A1
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG2A
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P11836
NCBI	P11836
Expiration Date	12 months from date of receipt.