


Product Datasheet

CXCL8 Biosimilar Antibody (orb1238042)

Description	CXCL8 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / IL-8 (Granulocyte chemotactic protein 1, Monocyte-derived neutrophil-activating peptide, IL-8 (1-77), IL-8 (9-77), Monocyte-derived neutrophil chemotactic factor, IL-8 (7-77), (Ser-IL-8)72, GCP/IL-8 protein V, MONAP, C-X-C motif chemokine 8, GCP/IL-8 protein I, Neutrophil-ac) [Homo sapiens]
Target	CXCL8
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P10145
NCBI	P10145
Expiration Date	12 months from date of receipt.



SDS-PAGE
Image