

Product Datasheet

ORF6 Antibody (orb1238085)

Description

ORF6 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

Anti-SARS-CoV-2 (COVID-19) ORF6 antibody (orb1238085) was raised against a peptide corresponding to 15 amino acids near the Carboxyl terminus of SARS-CoV-2 (COVID-19) ORF6 protein. The immunogen is located in the last 50 amino acids of SARS-CoV-2 (COVID-19) ORF6.

Target

ORF6

Preservatives

SARS-CoV-2 (COVID-19) ORF6 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) ORF6 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[P0DTC6](#)

NCBI

[QKU37157.1](#)

Dilution Range

Antibody validated: SARS-CoV-2 (COVID-19) ORF6 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA.

Expiration Date

12 months from date of receipt.