

Product Datasheet

Zika Envelope Antibody (orb1238086)

Description

Zika Envelope Antibody

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Zika Envelope antibody was raised against a peptide corresponding to 14 amino acids near the amino terminus of human Zika Envelope.
Target	Zika Envelope
Preservatives	Zika Envelope Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Zika Envelope antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Zika Envelope antibody can be used for Western blot at 1 - 2 µg/mL.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q32ZE1
NCBI	AMR39836
Dilution Range	Zika Envelope antibody can be used for Western blot at 1 - 2 µg/mL.
Expiration Date	12 months from date of receipt.