


Product Datasheet

ZMYM4 Antibody (orb1238131)

Description

ZMYM4 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	ZMYM4 antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of human ZMYM4. The immunogen is located within the last 50 amino acids of ZMYM4.
Target	ZMYM4
Preservatives	ZMYM4 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	ZMYM4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	ZMYM4 antibody can be used for detection of ZMYM4 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	167 kDa
Uniprot ID	Q5VZL5
NCBI	AAI27114
Dilution Range	ZMYM4 antibody can be used for detection of ZMYM4 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of ZMYM4 in human
t...