

Product Datasheet

TRPV1 Antibody (orb1238158)

Description

TRPV1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Rabbit polyclonal TRPV1 antibody was raised against a 16 amino acid peptide near the carboxy terminus of human TRPV1. The immunogen is located within the last 50 amino acids of TRPV1.

Target

TRPV1

Preservatives

TRPV1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TRPV1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

TRPV1 antibody can be used for detection of TRPV1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 92 kDa

Uniprot ID

[Q8NER1](#)

NCBI


[NP_542437](#)

Dilution Range

TRPV1 antibody can be used for detection of TRPV1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of TRPV1 in K562 c...